

CD3 (ABT-CD3) mouse mAb

Catalog No: #33204



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CD3 (ABT-CD3) mouse mAb
Brief Description	Mouse Monoclonal Ab
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	IHC;IF
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of human CD3. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Immunogen Description	Synthesized peptide derived from human CD3
Target Name	CD3E T3E
Uniprot	P07766
GeneID	916
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

IHC-p 1:100-500oO WB 1:500-2000. IF 1:50-200

Images



Western blot analysis of CD3Antibody at 1:1000 dilution.



Immunohistochemical analysis of paraffin-embedded Spleen-high magnification.

Note: This product is for in vitro research use only