

NPY1-R Polyclonal Antibody

Catalog No: #46415

Package Size: #46415-1 50ul #46415-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NPY1-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	NPY1-R Polyclonal Antibody detects endogenous levels of NPY1-R protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human NPY1R.
Target Name	NPY1R
Other Names	NPY1R; NPYR; NPY1; Neuropeptide Y receptor type 1; NPY1-R
Uniprot	P25929
GeneID	4886
Calculated MW	48kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000.

Images



Western Blot analysis of 3T3 cells using NPY1-R Polyclonal Antibody

Note: This product is for in vitro research use only