

# Goat Lactoferrin (LTF) ELISA Kit

Catalog No: #EK20512



Package Size: #EK20512-1 48T #EK20512-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Lactoferrin (LTF) ELISA Kit
Brief Description	ELISA Kit
Host Species	Goat
Applications	ELISA

## Application Details

Detect Range:0.781-50 ng/mL

Sensitivity:0.32ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only