

Rabbit Anti-MBOAT1 antibody

Catalog No: #47795



Package Size: #47795-1 50ul #47795-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-MBOAT1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human MBOAT1: 151-250/495
Target Name	MBOAT1
Other Names	1 acylglycerophosphoethanolamine O acyltransferase; LPE acyltransferase 1; LPEAT1; LysoPE acyltransferase; lysophosphatidylethanolamine acyltransferase 1; lysophospholipid acyltransferase; LysoPS acyltransferase; membrane bound O acyltransferase domain containing 1; O acyltransferase (membrane bound) domain containing 1; O acyltransferase domain containing protein 1; OACT1.
Uniprot	Q6ZNC8
GeneID	154141
Calculated MW	57kD
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:100-500 IHC-F=1:100-500 ICC=1:100-500 IF=1:100-500

Note: This product is for in vitro research use only