

SIRP- $\alpha$ 1 Polyclonal Antibody

Catalog No: #48134

Package Size: #48134-1 50ul #48134-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SIRP- $\alpha$ 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	SIRP- $\alpha$ 1 Polyclonal Antibody detects endogenous levels of SIRP- $\alpha$ 1 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human Sirp alpha1.
Target Name	SIRPA
Uniprot	P78324
GeneID	140885
Calculated MW	55kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000. IF 1:50-200

## Images



Western blot analysis of lysates from HepG2 cells, using Sirp alpha1 Antibody.

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using Sirp alpha1 Antibody.



---

Note: This product is for in vitro research use only