

CD69 Monoclonal Antibody

Catalog No: #48135

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD69 Monoclonal Antibody
Brief Description	Mouse Monoclonal Ab
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification
Applications	WB;IHC;IF;FCM;ELISA
Species Reactivity	Human
Specificity	CD69 Monoclonal Antibody detects endogenous levels of CD69 protein.
Immunogen Description	Purified recombinant fragment of human CD69 expressed in E. Coli.
Target Name	CD69
Uniprot	Q07108
GenelD	969
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Storage	Store at -20°C/1 year


Application Details

WB: 1/500 - 1/2000. IHC: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000.IF 1:50-200

Images



Western Blot analysis using CD69 Monoclonal Antibody against, Jurkat (1), L1210 (2) and TPH-1 (3) cell lysate.



Immunohistochemistry analysis of paraffin-embedded human Tonsil tissues with AEC staining using CD69 Monoclonal Antibody.

Note: This product is for in vitro research use only