

## GSDMD (N terminal) Rabbit pAb

Catalog No: #55651

Package Size: #55651-1 50ul #55651-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GSDMD (N terminal) Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF/ICC
Species Reactivity	Human
Crossing Reactivity	Human, Mouse, Rat
Immunogen Description	Recombinant fusion protein of human GSDMD(NP_079012.3).
Other Names	DF5L;DFNA5L;FKSG10;GSDMDC1;GSDMD
Uniprot	P57764
GeneID	79792
Calculated MW	53kDa
Formulation	Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB, 1:500 - 1:2000|IF/ICC, 1:50 - 1:200

## Images



Immunofluorescence analysis of A-431 cells using GSDMD (Full Length+N terminal) Rabbit pAb at dilution of 1:200 (40x lens). Blue: DAPI for nuclear staining.

## Background

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

---

Note: This product is for in vitro research use only