

# Human Carbonic anhydrase-related protein 11 (CA11) ELISA Kit



Catalog No: #EK20519

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20519-1 48T #EK20519-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Carbonic anhydrase-related protein 11 (CA11) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.78-50U/mL

Sensitivity:0.29U/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only