

Mouse Sex-determining region Y protein (SRY) ELISA Kit



Catalog No: #EK20525

Orders: order@signalwayantibody.com

Package Size: #EK20525-1 48T #EK20525-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Sex-determining region Y protein (SRY) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:14.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only