

## ZN141 rabbit pAb

Catalog No: #48287

Package Size: #48287-1 50ul #48287-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZN141 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human;Rat;Mouse
Specificity	This antibody detects endogenous levels of ZN141 at Human
Immunogen Description	Synthesized peptide derived from human ZN141
Target Name	ZN141
Uniprot	Q15928
GeneID	7700
Calculated MW	52kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1αO 500-2000

## Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only