

# Goat anti Rabbit/Mouse polymer HRP, Ready to Use

Catalog No: #BR009



Package Size: #BR009-1 5ml #BR009-2 10ml #BR009-3 100ml #BR009-4 1L

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat anti Rabbit/Mouse polymer HRP, Ready to Use
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC;IF
Species Reactivity	Mouse Rabbit
Specificity	This antibody detects Mouse or Rabbit IgG antibodies with labeled HRP
Storage	Store at 4°C/1 year

## Application Details

Ready to Use

Note: This product is for in vitro research use only