

FMR1 Monoclonal Antibody

Catalog No: #48292

Package Size: #48292-1 50ul #48292-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FMR1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human
Specificity	FMR1 Monoclonal Antibody detects endogenous levels of FMR1 protein.
Immunogen Description	Purified recombinant fragment of human FMR1 expressed in E. Coli.
Target Name	FMR1
Other Names	FMR1; Fragile X mental retardation protein 1; FMRP; Protein FMR-1
Uniprot	Q06787
GeneID	2332
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000.


Images



Western Blot analysis using FMR1 Monoclonal Antibody against Jurkat (1) and K562 (2) cell lysate.



Immunofluorescence analysis of NIH/3T3 cells using FMR1 Monoclonal Antibody (green).



Immunohistochemistry analysis of paraffin-embedded human brain tissues

Note: This product is for in vitro research use only