

# Recombinant Human Delta-type opioid receptor(OPRD1)

Catalog No: #AP89656

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89656-1 20ug

## Description

Product Name	Recombinant Human Delta-type opioid receptor(OPRD1)
Brief Description	Recombinant Proteins
Accession No.	P41143
Uniprot	P41143
Target Species	Homo sapiens (Human)
Calculated MW	41.9 kDa
Tag Info	N-terminal 10xHis-tagged
Target Sequence	Protein Length:Full Length Expression Region: 1-372aa MEPAPSAGAE LQPPLFANASDAYPSACPSAGANASGPPGARSASSLALAI AITALYSAVCAVGLLG NVLVMFGI VRYTKMKTATNIYIFNLALADALATSTL PFQSAKYL METWPF GELLCKAVLSIDYYNMFTSIFTLTMM SVDRYIAV CHPVKALDFRTPAKAKLINICIWV LASGVGVPI MVMAVTRPRDGA VVCM LQFPSPSWYWDTVTKICVFLFAFV PILIITVCYGLMLLR LRSVRLLSGSKEKDRSLRRITRMVLVVGFVVCWAPIHIFVIVWTLVDIDRRDPLVVAALH LCIALGYANSSLN PVLYAFLDENFKRCFRQLCRKPCGRPDSSFSRAREATARERTACTPSDGPGGGAAA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only