

Anti-WDR61 Antibody

Catalog No: #48328

Package Size: #48328-1 50ul #48328-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-WDR61 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB;IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human WDR61.
Target Name	WDR61
Other Names	WD repeat-containing protein 61; Meiotic recombination REC14 protein homolog; SKI8 homolog; Ski8
Accession No.	Q9GZS3
Uniprot	Q9GZS3
GeneID	80349
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of WDR61 expression in HeLa (A), mouse brain (B), mouse kidney (C), rat pancreas (D) whole cell lysates.



Immunohistochemical analysis of WDR61 staining in human placenta formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Note: This product is for in vitro research use only