

Monkey Immunoglobulin G (IgG) ELISA Kit

Catalog No: #EK20529



Package Size: #EK20529-1 48T #EK20529-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monkey Immunoglobulin G (IgG) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range: 3.12-200ug/ml

Sensitivity: 1.36ug/ml

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: competitive inhibition

Note: This product is for in vitro research use only