

NPRL3 rabbit pAb

Catalog No: #48618

Package Size: #48618-1 50ul #48618-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NPRL3 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse
Specificity	This antibody detects endogenous levels of NPRL3 at Human/Mouse
Immunogen Description	Synthesized peptide derived from human NPRL3
Target Name	NPRL3
Other Names	Nitrogen permease regulator 3-like protein (-14 gene protein) (Alpha-globin regulatory element-containing gene protein) (Protein CGTHBA)
Uniprot	Q12980
GeneID	8131
Calculated MW	60kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.170% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only