

Mouse Interleukin 1 Receptor Associated Kinase 4 (IRAK4) ELISA Kit

Catalog No: #EK20542

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20542-1 48T #EK20542-2 96T

Description

Product Name	Mouse Interleukin 1 Receptor Associated Kinase 4 (IRAK4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000 pg/mL

Sensitivity:6.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only