

Mouse Glutathione Peroxidase 1 (GPX1) ELISA Kit

Catalog No: #EK20558



Package Size: #EK20558-1 48T #EK20558-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Glutathione Peroxidase 1 (GPX1) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Mouse |

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.138ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only