

# Mouse Carnitine palmitoyltransferase 2, mitochondrial (CPT2) ELISA Kit



Catalog No: #EK20545

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20545-1 48T #EK20545-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Carnitine palmitoyltransferase 2, mitochondrial (CPT2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:62.5-4000 pg/mL
Sensitivity:24.5pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only