

# Mouse Carnitine Acylcarnitine Translocase (CACT) ELISA Kit



Catalog No: #EK20559

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20559-1 48T #EK20559-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Carnitine Acylcarnitine Translocase (CACT) ELISA Kit
Brief Description	ELISA Kit

## Application Details

Detect Range: 0.156-10 ng/mL

Sensitivity: 0.067 ng/ml

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only