

## GPR84 Polyclonal Antibody APC Conjugated

Catalog No: #C90969APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR84 Polyclonal Antibody APC Conjugated
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt,IF,ICC
Species Reactivity	Human,Mouse
Immunogen Description	KLH conjugated synthetic peptide derived from human GPCR EX33 protein
Target Name	GPR84
Other Names	G protein coupled receptor 84; G-protein coupled receptor 84; G protein coupled receptor EX33; GPCR 4; GPCR4; GPR 84; GPR84_HUMAN; GPR84 (EX33); EX33; Inflammation related G protein coupled receptor EX33; Inflammation-related G-protein coupled receptor EX33; GPCR EX33; EX 33; GPCR GPR84;
Uniprot	Q9NQS5
GeneID	53831
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt 1&amp;mu;g/Test.IF 1:100-500.ICC 1:100-500

Note: This product is for in vitro research use only