

NRIP1 Rabbit pAb

Catalog No: #55653

Package Size: #55653-1 50ul #55653-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NRIP1 Rabbit pAb
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Recombinant fusion protein of human NRIP1(NP_003480.2).
Other Names	NRIP1;RIP140
Uniprot	P48552
GeneID	8204
Target Species	Human
Calculated MW	126kDa
SDS-PAGE MW	150KDa
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Application Details

WB, 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse brain, using NRIP1 Polyclonal Antibody at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit.

Exposure time: 3min.

Background

Nuclear receptor interacting protein 1 (NRIP1) is a nuclear protein that specifically interacts with the hormone-dependent activation domain AF2 of nuclear receptors. Also known as RIP140, this protein modulates transcriptional activity of the estrogen receptor.

Note: This product is for in vitro research use only