

# Streptavidin Cy5 Conjugated

Catalog No: #LS1014



Package Size: #LS1014-1 100ul #LS1014-2 1ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Streptavidin Cy5 Conjugated
Brief Description	Streptavidin Cy5 Conjugated
Purification	AC
Conjugates	Cy5
Other Names	SA; Streptavidin / Cy5
Concentration	σO 0.5 mg/ml
Formulation	10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only