

Human BTN3A2/ Butyrophilin subfamily 3 member A2 ELISA Kit



Catalog No: #EK21111

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human BTN3A2/ Butyrophilin subfamily 3 member A2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Uniprot	P78410

Application Details

Detect Range:78-5000pg/mL

Sensitivity:28pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only