

Mouse Pdgfra/ Platelet-derived growth factor receptor alpha ELISA Kit



Catalog No: #EK21124

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Pdgfra/ Platelet-derived growth factor receptor alpha ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Mouse |
| Uniprot | P26618 |

Application Details

| |
|--|
| Detect Range:0.312-20ng/mL |
| Sensitivity:0.171ng/mL |
| Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Volume: 50 μL |
| Assay Time:1-4.5h |
| Detection wavelength:450 nm |

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only