

IL18 binding protein Antibody

Catalog No: #57170

Package Size: #57170-1 50ul #57170-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	IL18 binding protein Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human,Rat
Specificity	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Immunogen Description	A synthesized peptide derived from human IL18 binding protein
Other Names	IL18 binding protein; IL18 Bpa; IL18BP; IL18BPa; Interleukin 18 binding protein; MC51L 53L 54L homolog gene product; Tadekinig alfa;
Uniprot	O95998
Calculated MW	40kDa
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500-1:2000

IHC 1:50-1:200

Images



Western blot analysis of IL18 binding protein in HeLa cell lysate.

Background

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

Note: This product is for in vitro research use only