

G6PC rabbit pAb

Catalog No: #48660



Package Size: #48660-1 50ul #48660-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G6PC rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB;IHC
Species Reactivity	Human, Mouse,Rat
Specificity	This antibody detects endogenous levels of G6PC at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human C-terminal G6PC
Target Name	G6PC G6PT
Other Names	Glucose-6-phosphatase (G-6-Pase) (G6Pase) (EC 3.1.3.9) (Glucose-6-phosphatase alpha) (G6Pase-alpha)
Uniprot	P35575
GeneID	2538
Calculated MW	39kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 IHC 1:50-200

Note: This product is for in vitro research use only