

GPR91 Polyclonal Antibody

Catalog No: #48677

Package Size: #48677-1 50ul #48677-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GPR91 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse;Rat;Pig
Specificity	GPR91 Polyclonal Antibody detects endogenous levels of GPR91 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human SUCNR1.
Target Name	SUCNR1
Other Names	SUCNR1; GPR91; Succinate receptor 1; G-protein coupled receptor 91; P2Y purinoceptor 1-like
Uniprot	Q9BXA5
GeneID	56670
Calculated MW	38kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. IF 1:200 - 1:1000. ELISA: 1:10000.

Images



Western blot analysis of lysates from HUVEC and MCF-7 cells, using SUCNR1 Antibody. The lane on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only