

## PIEZ1 rabbit pAb

Catalog No: #48683

Package Size: #48683-1 50ul #48683-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PIEZ1 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC;WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of PIEZ1 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human C-terminal PIEZ1
Target Name	PIEZO1 FAM38A KIAA0233
Other Names	Piezo-type mechanosensitive ion channel component 1 (Membrane protein induced by beta-amyloid treatment) (Mib) (Protein FAM38A)
Uniprot	Q92508
GeneID	9780
Calculated MW	277kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000 IHC 1:50-200

Note: This product is for in vitro research use only