

MouseB Lipoprotein, B aB (Lpa)B ELISAB Kit

Catalog No: #EK20555



Package Size: #EK20555-1 48T #EK20555-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MouseB Lipoprotein, B aB (Lpa)B ELISAB Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range: 0.6ng/ml - 7.2ng/ml

Sensitivity: 0.12ng/ml

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only