

# Fish Interleukin 1 Beta (IL1b) ELISA Kit

Catalog No: #EK20556



Package Size: #EK20556-1 48T #EK20556-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fish Interleukin 1 Beta (IL1b) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Fish

## Application Details

Detect Range: 6pg/ml-72pg/ml

Sensitivity: 1.2pg/ml

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only