

# Fish Tumor Necrosis Factor Alpha(TNF $\alpha$ )ELISA Kit

Catalog No: #EK20560



Package Size: #EK20560-1 48T #EK20560-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fish Tumor Necrosis Factor Alpha(TNF $\alpha$ )ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Fish

## Application Details

Detect Range:7.5pg/ml-90pg/ml

Sensitivity:1.2pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only