

Recombinant Human Pachytene checkpoint protein 2 homolog(TRIP13)

Catalog No: #AP89664

Orders: order@signalwayantibody.com

Package Size: #AP89664-1 20ug #AP89664-2 100ug #AP89664-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Pachytene checkpoint protein 2 homolog(TRIP13)
Brief Description	Recombinant Protein
Target Name	His-tagged
Uniprot	Q15645
Target Species	Homo sapiens (Human)
Calculated MW	54.6 kDa
Tag Info	N-terminal 10xHis-tagged
Target Sequence	<p>Protein Length:Full Length</p> <p>Expression Region:1-432aa</p> <p>MDEAVGDLKQALPCVAESPTVHVEVHQSGSSTAKKEDINLSVRKLLNRHNIVFGDYTWTEFDEPFLTRNVQSV SIIDTELKVKDSQPIDLSACTVALHIFQLNEDGPSENLEEETENIIAANHWWLPAAEFHGLWDSLVDVEVKSH LLDYVMTLLFSDKNVNSNLITWNRVLLHGPPGTGKTSLCKALAQKL TIRLSSRYRYGQLIEINSHSLFSKWFS ESGKLVTKMFQKIQLDIDDKDALVFLIDEVESLTAARNACRAGTEPSDAIRVVNAVLTQIDQIKRHSNVILTTS NITEKIDVAFVDRADIKQYIGPPSAAAIFKIYLSCLEELMKCQIYPRQQLTLRELEMIGFIENNVSKLSLLLNDISR KSEGLSGRVLRKLPFLAHALYVQAPTVTIEGFLQALSLAVDKQFEERKKLAAYI</p>
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only