

MRC1 Antibody APC-Cy7 Conjugated

Catalog No: #C05371APC-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MRC1 Antibody APC-Cy7 Conjugated
Brief Description	rabbit polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human MRC1/CD206
Conjugates	APC-Cy7
Target Name	MRC1
Other Names	Mannose Receptor; C-type lectin domain family 13 member D; CD 206; CD206; CD206 antigen; CLEC13D; CLEC13DL; Macrophage mannose receptor; Mannose receptor C type 1; MMR; MRC 1; MRC-1; MRC1_HUMAN; Macrophage mannose receptor 1; MRC1L1; C-type lectin domain family 13 member D-like; Macrophage mannose receptor 1-like protein 1.
Uniprot	P22897
GeneID	4360
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt:2&mu;g /Test

Note: This product is for in vitro research use only