

Rat Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit

Catalog No: #EK20564



Package Size: #EK20564-1 48T #EK20564-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Carbohydrate Antigen 15-3 (CA15-3) ELISA Kit
Brief Description	ELISA Kit
Host Species	Rat
Purification	ELISA

Application Details

Detect Range:7.81-500 mIU/mL

Sensitivity:3.2mIU/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only