

Canine Fibroblast Growth Factor 23 (FGF23) ELISA Kit



Catalog No: #EK20568

Orders: order@signalwayantibody.com

Package Size: #EK20568-1 48T #EK20568-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Canine Fibroblast Growth Factor 23 (FGF23) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range:7.81-500 pg/ml

Sensitivity:3.6pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only