

# Mouse Hepatocyte cell adhesion molecule (hepacam) ELISA Kit



Catalog No: #EK20569

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20569-1 48T #EK20569-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Hepatocyte cell adhesion molecule (hepacam) ELISA Kit
Brief Description	ELISA Kit
Host Species	Mouse
Purification	ELISA

## Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:27pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only