

Mouse Ovalbumin IgG (OVA-IgG) ELISA Kit

Catalog No: #EK20573



Package Size: #EK20573-1 48T #EK20573-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ovalbumin IgG (OVA-IgG) ELISA Kit
Host Species	Mouse
Purification	ELISA

Application Details

Detect Range:31.2-2000 pg/mL

Sensitivity:12.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only