

Porcine Alpha-1-Microglobulin/Bikunin Precursor (a1M) ELISA Kit



Catalog No: #EK20574

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20574-1 48T #EK20574-2 96T

Description

Product Name	Porcine Alpha-1-Microglobulin/Bikunin Precursor (a1M) ELISA Kit
Brief Description	ELISA Kit
Purification	ELISA
Species Reactivity	Pig

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.37ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only