

# Mouse Ifnar1/ Interferon alpha/beta receptor 1 ELISA Kit



Catalog No: #EK21125

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK21125-1 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ifnar1/ Interferon alpha/beta receptor 1 ELISA Kit
Brief Description	ELISA Kit
Host Species	Mouse
Purification	ELISA
Other Names	Type I interferon receptor 1, IFN-R-1, Ifar, Ifnar

## Application Details

Detect Range:78-5000pg/mL

Sensitivity:39pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only