

Mouse Ifnar2/ Interferon alpha/beta receptor 2 ELISA Kit



Catalog No: #EK21126

Orders: order@signalwayantibody.com

Package Size: #EK21126-1 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ifnar2/ Interferon alpha/beta receptor 2 ELISA Kit
Brief Description	ELISA Kit
Host Species	Mouse
Purification	ELISA
Other Names	Type I interferon receptor 2, IFN-R-2

Application Details

Detect Range:78-5000pg/mL

Sensitivity:40.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only