

Mouse Ifngr1/ Interferon gamma receptor 1 ELISA Kit

Catalog No: #EK21127



Package Size: #EK21127-1 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ifngr1/ Interferon gamma receptor 1 ELISA Kit
Brief Description	ELISA Kit
Host Species	Mouse
Purification	ELISA
Other Names	Interferon gamma receptor alpha-chain, IFN-gamma-R-alpha, CD119, IFN-gamma receptor 1, Ifngr

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.156ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only