

# Mouse serum albumin Polyclonal Antibody HRP Conjugated

Catalog No: #C91594HRP

Package Size: #C91594HRP-1 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse serum albumin Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	IHC ELISA
Species Reactivity	Mouse,Human
Immunogen Description	Mouse serum albumin purified from mouse serum:full length
Conjugates	HRP
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC 1:100-500 ELISA 1:5000-10000

Note: This product is for in vitro research use only