

CD19 Monoclonal Antibody

Catalog No: #48814

Package Size: #48814-1 50ul #48814-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD19 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human
Specificity	CD19 Monoclonal Antibody detects endogenous levels of CD19 protein.
Immunogen Description	Purified recombinant fragment of human CD19 expressed in E. Coli.
Target Name	CD19
Other Names	CD19; B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; Differentiation antigen CD19; T-cell surface antigen Leu-12; CD antigen CD19
Uniprot	P15391
GeneID	930
Concentration	1 mg/ml
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. IHC 1:200 - 1:1000. ELISA: 1:10000.. IF 1:50-200

Images



Immunohistochemistry analysis of paraffin-embedded human normal lymph node, showing cytoplasmic localization with DAB staining using CD19 Monoclonal Antibody.

Note: This product is for in vitro research use only