

# Human Cyclic AMP-responsive element-binding protein 3-like protein 4 $\alpha$ O CREB3L4 $\alpha$ O ELISA Kit

Catalog No: #EK20586

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20586-1 48T #EK20586-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Cyclic AMP-responsive element-binding protein 3-like protein 4 $\alpha$ O CREB3L4 $\alpha$ O ELISA Kit
Brief Description	ELISA Kit
Purification	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.156-10ng/ml

Sensitivity:0.057ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only