

## MSP Rabbit mAb

Catalog No: #57200

Package Size: #57200-1 50ul #57200-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MSP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Serine/threonine-protein kinase 4
Other Names	serine threonine kinase 3/4; KRS1; MST2/KRS2; MST1; YSK3; TIAC
Accession No.	Swiss-Prot:Q13043 Gene ID:6789
Uniprot	Q13043
GeneID	6789;
Calculated MW	56 kDa
SDS-PAGE MW	60 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details


IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

## Images



Western blot analysis of MSP/MST1 in K562, rat Brain, C6, 3T3, HeLa lysates using MSP/MST1 antibody.

Immunocytochemistry analysis of MSP (green) in K562 using MSP antibody, and DAPI(blue).



---

Note: This product is for in vitro research use only