

PTF1A Rabbit mAb

Catalog No: #57203

Package Size: #57203-1 50ul #57203-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PTF1A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of mouse PTF1A
Other Names	Pancreas-specific transcription factor 1a; bHLH transcription factor p48; PTF1-p48
Accession No.	Swiss-Prot:Q9QX98 Gene ID:19213
Uniprot	Q9QX98
GeneID	19213;
Calculated MW	35 kDa
SDS-PAGE MW	48 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of PTF1A in K562, C6 lysates using PTF1A antibody.

Note: This product is for in vitro research use only