

MAGEC2 Rabbit mAb

Catalog No: #57204

Package Size: #57204-1 50ul #57204-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MAGEC2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human MAGEE1
Other Names	CT10
Accession No.	Swiss-Prot:Q9UBF1 Gene ID:51438
Uniprot	Q9UBF1
GeneID	51438;
Calculated MW	41 kDa
SDS-PAGE MW	55 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of MAGEE1 in K562 lysates using MAGEE1 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse testis using MAGEC2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note: This product is for in vitro research use only