

PBR Rabbit mAb

Catalog No: #57205



Package Size: #57205-1 50ul #57205-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PBR Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human PBR
Other Names	TSPO; BZRP; MBR; Translocator protein; Mitochondrial benzodiazepine receptor; PKBS; Peripheral-type benzodiazepine receptor; PBR
Accession No.	Swiss-Prot:B1AH88 Gene ID:706
Uniprot	B1AH88
GeneID	706;
Calculated MW	11 kDa
SDS-PAGE MW	19 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200


Images



Western blot analysis of PBR/TSPO in 3T3, Hela lysates using PBR antibody



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PBR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PBR/TSPO (green) in Hela using PBR/TSPO antibody, and DAPI (blue)

Note: This product is for in vitro research use only