

## TGF beta 1 Rabbit mAb

Catalog No: #57206

Package Size: #57206-1 50ul #57206-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TGF beta 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of mouse TGF beta 1
Other Names	TGF beta 1; TGFB; CED; LAP
Accession No.	Swiss-Prot:P04202 Gene ID:21803
Uniprot	P04202
GeneID	21803;
Calculated MW	44 kDa
SDS-PAGE MW	44 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/500-1/1000

## Images



Western blot analysis of TGF beta 1 in Raw264.7 lysates using TGF beta 1 antibody.



Western blot analysis of TGF beta 1 in Jurkat, C6 lysates using TGF beta 1 antibody.

---

Note: This product is for in vitro research use only